



YINGCAI BEAUTY

PROFESSIONAL BEAUTY EQUIPMENT FACTORY

Patented Products. Global Partners



PLUS ONE

TA405

JUST FOR PROFESSIONAL



Laser Hair Removal

755nm/808nm/1064nm/940nm

- STRONG MAGNETIC HEADS
- 6 SPOT SIZE FOR CHOICE
- FOUR WAVELENGTHS




ENERGY LOSS COMPARISON



6 SPOT SIZE FOR CHOICE




DIFFERENT SPOTS CORRESPOND TO DIFFERENT PARTS




AUTOMATIC SPOT RECOGNITION




CORE CONFIGURATION



Powerful configuration enables the machine to work continuously for 24 hours



FOUR WAVELENGTHS FOR ALLTYPE SKIN HAIR REMOVAL



5000W

STRONG POWER

4 Handles


Muscle Building

EMS 15 Tesla

- Easily thin body
- Get into shape quickly
- No movement



Accessory



Universal Handle*2
For ABS,Thighs,Belly and Buttocks



Limb handles*2
For Arms,Lights,Thighs and Shank

Optional



Therapeutic cushions
For postpartum recovery, pelvic floor muscles training and private treatment



Customizable treatment plans

Can be tailored to each individual's specific goals and body type. The treatment can target various areas of the body, including the arms, abdomen, buttocks, thighs and calves.

Abdominal handle



Arm handle



SEAT-TYPE WORK HANDLE (OPTIONAL)



Pelvic floor training, postpartum recovery, reproductive privacy.
Hip lifting effect

Configuration

	Standard	Upgrade	Luxury
Control Panel			
Capacitor			
	One Controls, Two Capacitors, Four Handles.	Two Controls, 4 Capacitors, Four Handles.	Four Controls, 4 Capacitors, Four Handles.

HIGH INTENSITY FOCUSED ELECTRO-MAGNETIC TECHNOLOGY

-50%



Can increase the efficiency of fat reduction

+40%

Can help promote better muscle growth



Amazing Effect



11 IN 1 Slimming Machine

BG16



Focused
ultrasound



Focused ultrasound(S)
Promote Metabolism And Soften Tissues



Face RF
Massage the face to improve wrinkles



Vacuum RF
Body Shaping, Cupping



Body Lift
Massage The Skin To Promote Metabolism



Cift Roller
Massage, Sculpt And Relax Muscles



Focused ultrasound(L)
Promote The Metabolism Of Excess Substances In The Body



Body RF
Firm And Shape The Abdomen





Pop-Up
Breast And Butt Lifting, Massage



Face Lift
Microcurrents Enhance Facial Elasticity



Microelectric



DPL PRO


Add three additional ANCE/VASCULAR/CALIBRATION



ELEVEN

Wavelength Filters Optional

- Legs
- Back and chest
- Intimates/ bikini line
- Face (above lips, nose, cheeks, chin, ears, eyebrows)
- Arms and underarms




11 Wavelength Filters

1 420nm <ul style="list-style-type: none"> ● Acne treatment ● Oil-control effect ● Anti inflammatory and sedative 	2 500nm <ul style="list-style-type: none"> ● Acne treatment ● Anti-inflammatory ● Whitening and improving skin tone 	3 515nm <ul style="list-style-type: none"> ● Epidermis pigmentation
4 550nm <ul style="list-style-type: none"> ● Pigmentation ● Uneven skin tone ● Minimize pores 	5 560nm <ul style="list-style-type: none"> ● Pigmentation ● Vascularisation and skin care 	6 590nm <ul style="list-style-type: none"> ● Dermis larger vessels ● Skir care ● Pigmentation and vascularization
7 615nm <ul style="list-style-type: none"> ● Larger vessels in the face ● General IPL skin care ● Pigmentation and vascularization 	8 640nm <ul style="list-style-type: none"> ● Anti-aging ● Deep repair ● Whitening and brightening the skin 	9 695nm <ul style="list-style-type: none"> ● Hair removal
10 ACNE 	11 VASCULAR 	

Precise Treatment

Pulse number, pulse width, and pulse interval time are adjustable



DIODE LASER+PICO LASER


1064+755+808nm

Hair Removal Handle :
Standard Configuration:


1200W

Optional Configurations:

600W 800W 1000W



Remove
Vigorous Hair




Freckle
removing



Improves red
vascular lesions



Enhances skin
appearance



Remove
tattoo

DIODE LASER



High Power

600W/800W/1000W/1200W



Handle Shots

100 million shots



COOLING SYSTEM

LOWEST TEMPERATURE
-42°C



PICO LASER

SUITABLE FOR ALL SKIN TYPES AND
REMOVE ALL COLORS TATTOO AND PIGMENTS



PICOSECOND

4 WAVELENGTHS

TATTOO REMOVAL MACHINE



2500
High Power



5 Treatment Probes



Remove Tattoos


Fine Pores

Skin Rejuvenation

Whitening Skin

MORE PROFESSIONAL MORE PRECISE

5 large probes for precise targeting of pigments



PICOSECOND HANDLE



2500W

Ensure stable performance



5 Million Shoots

More than 5 million shoots used long time



1500mj

High energy ensure the obvious treatment effects



Optional probe spot

Choice of honeycomb focusing probe or dot focusing probe



4 Wavelengths

Suitable for all skin types Remove all tattoos and pigment



HZ

10HZ

10Hz ensure the treatment more effective

PERMANENT HAIR REMOVAL

2 IN 1

Insulated Electric Needles
+ Hot And Cold Hammers



Hot And Cold Hammers



Hair Removal



Hair Removal



Treat The Skin



Photo Rejuvenation

HANDLE DISPLAY



Hot and cold hammer


After completing hair removal, use an ice hammer to cool down the hair removal area



Insulated Electric Needles

The effect of removing hair is achieved by breaking hair follicles through the action of electric current.


TECHNICAL EXPLANATION



Prone to scarring

Traditional electroacupuncture technique:

The entire needle body is conductive, so while it damages hair follicles, it also causes **damage** to the epidermis and is **prone to scarring**.



Without damaging the epidermis

Insulated electroacupuncture technology:

Insulated electroacupuncture technology: using specially designed insulated needles, **only the needle tip part conducts electricity**. it effectively destroys hair follicles **without damaging the epidermis**.

History of Electrolysis Hair Removal



Version 1.0 2.0MHZ



Version 2.0 13.5MHZ




Version 3.0 27.5MHZ



Version 4.0 36.5MHZ

5 IN 1 Permanent

Multifunctional skin problem treatment



Diameter 1 mm

Other

x20 boxes A box of 50 needles

Diameter 0.02 mm

Advantage

- 1. Medical-grade
- 2. 24k gold plated
- 3. Achieved insulation that is non-invasive and does not harm the skin

A whole box contains 1000 needles
K2 / K3 / K4 / K5

Hot and cold hammer

After completing hair removal, use an ice hammer to cool down the hair removal area

Insulated Electric Needles

The effect of removing hair is achieved by breaking hair follicles

Plasma pen

With specially designed needles conducted only a part of the energy, remove spots, pigmentation, freckles, etc.

Hair removal with tweezers

With specially designed tweezers conducted only a part of the electric

Anti-static wrist strap

Wearing it on your hand during operation can effectively relieve pain.

7D HIFU + RF Microneedle + Vaginal

HB08



Germany Ultrasonic
film Germany



Strong energy,
more significant effect



1.5mm 2.0mm 3.0mm 4.5mm



6.0mm 9.0mm 13.0mm



64Pins

64 Pins apply to body



25 Pins

25 Pins apply to face



10Pins

10 Pins suitable for
dead corner



Nano


Nano chip apply to face



3.0mm



4.5mm



Diode Laser+Cooling

New

- 1 New appearance design
- 2 Remote control system
- 3 USA coherent laser bar
- 4 High power TEC condenser



Target Imported From Germany

755nm/808nm/940nm/1064nm




100 Million Shots

Better lasers can treat more patients

600W 800W 1000W 1200W 1600W

Newest Technology Fast Axis Collimator

Exclusively customized coated crystal, output light energy is more concentrated: equipped with laser calibration to make the energy more concentrated




6 Different Big Spot Size

Magnetic striker embedded design, equipped with dust cover, safer and more durable



SIZE:



Screen linkage Improve working efficiency



CO2

Fractional Laser Machine



Patented CO2 Laser

Function

 Skin Renewing And Resurfacing	 Skin Tightening Wrinkle Removal
 Pore Minimization	 Sun Damage Recovery
 Acne Removal	 Remove The Pigment
 Smooth Scars And Dermabrasion	 Vaginal Tightening Vulve Rejuvenatio


Advantages

- Laser Power : **60w**
- Display : **15.6-Inch**
- Light Guide Arm : **6 Joints**-Korea Imports
- Cooling System : Water Cooling
- Wavelength : **10600nm**
- Mode : **4 Working Modes**
(Fractional mode, Pulse mode, Gynae mode, Vulva mode)
- Focus Spot Diameter : 20*20mm

Machine Functions

Variety Of Scan Shapes Available

1 Fractional Co2 Laser Probe
For Skin Rejuvenation



Treatment Graphic

Suitable For Many Different Areas Of Treatment.



Treatment Head




1 Normal Co2 Laser Probe
For Surgery Cutting

- ✓ Safe Precision
- ✓ No Downtime Comfortable

3 Vaginal probes

For vaginal skin rejuvenation

- A.** Vulva Rejuvenation Probe
- B.** Inner Vagina Rejuvenation Probe
- C.** Treatment of Cervical Erosion Probe



 Skin renewing and resurfacing	 Skin tightening wrinkle removal	 Smooth scars and dermabrasion	 Sun damage recovery
 Acne removal	 Remove the pigment	 Pore minimization	 Vaginal tightening, Vulve Rejuvenation.



Global Expertise, Focused on Your Challenges.

With a proven track record in over 50 countries, Yingcai has consistently delivered where it matters most. We have successfully executed complex projects for leading clients across the **United States, Germany, France, and the United Kingdom**, tackling their most demanding aesthetic technology challenges. This depth of experience means we are uniquely equipped to ensure your project's success, no matter the complexity. When you work with Yingcai, you tap into a world of proven solutions.



YINGCAI BEAUTY

PROFESSIONAL BEAUTY EQUIPMENT FACTORY

Company Name:Guangzhou YingCai Beauty Machine Co.,Ltd

Contact Person:Spring

WhatsAPP/Wechat:+86 13250766080

Add:6F# No.5 Yingbian West Road,Baiyun District,Guangzhou

Guangdong, 510000,China

<http://ycblaser.com>